

Intercept® Oral Fluid Drug Testing



FDA-Cleared, Laboratory-Based Oral Fluid Drug Testing.

Lower Costs with Increased Reliability




Maximize Productivity and Hiring Time

- Hire employees faster
- Keep employees on the job, minimize lost work time
- Avoid scheduling inconveniences and fees of urine collection
- Eliminate the potential time wasted by shy bladder syndrome and the necessity to maintain a specialized restroom

Deliver Oral Fluid Convenience

- Speeds up and controls collection process
- Streamline collections done anytime, anywhere
- Preserves employee dignity with non-invasive collection process
- Eliminates adulteration and dilution tactics with directly observed collections every time

Intercept® Speeds Up the Drug Testing Process



Intercept® – Laboratory Confidence With Oral Fluid Convenience

	INTERCEPT®	Urine
Days to results for applicant testing ¹	1 day	2-5 days
Time off the job for collection (average) ²	Under 15 minutes	30-90 minutes
Ecstasy confirmation included in standard panel (NIDA-5) ³	Yes	No
Prevents tampering, dilution, or masking agents ⁴	Yes	No
Certified laboratory analysis	Yes	Yes
Training and certification for collection personnel ⁵	<1 hour, CD-based, very easy for existing employees	1-day seminar, many variables on diluted, adulterated or substituted samples

Swab



1. Place pad between cheek and gums for 5 minutes.

Snap



2. Insert pad into open vial, then snap wand.

Seal



3. Recap the vial and apply seal. Send to laboratory with chain of custody form.

For more information call 800.672.7873, or visit us online: ChooseIntercept.com

The Value of a Reliable and Trusted Brand

Intercept®

Oral Fluid Drug Test

1. Pre-employment urine drug screen collections and laboratory result reporting times will vary according to individual corporate policies and laboratory arrangements.
2. Oral Fluid Symposium: Advancements and Applications of Oral Fluid Drug Testing; Addiction Professional; Sept, 2004.
3. OraSure Amphetamine-Specific Intercept® MICRO-PLATE EIA Package Insert.
4. Oral fluid drug tests: Effects of adulterants and foodstuffs; Wong, Tran, Tung; Forensic Science International; Volume 150, Issue 2, 10, June 2005; 175-180.
5. Urine collection training times may vary by laboratory vendor.



OraSure Technologies, Inc.

220 East First Street, Bethlehem, PA 18015 USA, phone: 800.ORASURE
web: www.OraSure.com • www.ChooseIntercept.com

©2008, 2009 OraSure Technologies, Inc. • Intercept® is a registered trademark of OraSure Technologies, Inc.
Item# INTO109 (Rev. 11/10)

